



dinB Antibody

Product Code	CSB-PA676401ZA01CBAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q487H6
Immunogen	Recombinant Colwellia psychrerythraea dinB protein
Raised In	Rabbit
Species Reactivity	Colwellia psychrerythraea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase IV (Pol IV) (EC 2.7.7.7), dinB
Immunogen Species	Colwellia psychrerythraea (strain 34H / ATCC BAA-681?(Vibrio psychroerythus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dinB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.