



GFP Antibody

Product Code	CSB-PA337004ZA01ADO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42212
Immunogen	Recombinant Aequorea victoria Green fluorescent protein (GFP) (1-238AA)
Raised In	Rabbit
Species Reactivity	Aequorea victoria
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Green fluorescent protein, GFP
Immunogen Species	Aequorea victoria (Jellyfish)
Lead Time	Made-to-order
Target Names	GFP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.