



hamp Antibody

Product Code	CSB-PA010124ZA01LLH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1Z0M0
Immunogen	Recombinant Larimichthys crocea hamp protein
Raised In	Rabbit
Species Reactivity	Larimichthys crocea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepcidin, hamp
Immunogen Species	Larimichthys crocea (Large yellow croaker) (Pseudosciaena crocea)
Lead Time	Made-to-order (12-14 weeks)
Target Names	hamp
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.