



Non-specific lipid-transfer Antibody

Product Code	CSB-PA020856ZA01AOH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C088
Immunogen	Recombinant Artemisia vulgaris Non-specific lipid-transfer protein
Raised In	Rabbit
Species Reactivity	Artemisia vulgaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-specific lipid-transfer protein (LTP) (Pollen allergen Art v 3) (allergen Art v 3) (Fragment)
Immunogen Species	Artemisia vulgaris (Mugwort)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.