



HN Antibody

Product Code	CSB-PA336075ZA01NCW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35741
Immunogen	Recombinant Newcastle disease virus HN protein
Raised In	Rabbit
Species Reactivity	Newcastle disease virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemagglutinin-neuraminidase (EC 3.2.1.18), HN
Immunogen Species	Newcastle disease virus (strain Her/33) (NDV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HN
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.