



ctxB Antibody

Product Code	CSB-PA360704ZA01VEX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01556
Immunogen	Recombinant Vibrio cholerae serotype O1 ctxB protein
Raised In	Rabbit
Species Reactivity	Vibrio cholerae serotype O1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholera enterotoxin subunit B (Cholera enterotoxin B chain) (Cholera enterotoxin gamma chain) (Cholera toxin B subunit), ctxB, toxB
Immunogen Species	Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ctxB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.