



CTAG1A Antibody

Product Code	CSB-PA006125ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78358
Immunogen	Recombinant Homo sapiens CTAG1A protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cancer/testis antigen 1 (Autoimmunogenic cancer/testis antigen NY-ESO-1) (Cancer/testis antigen 6.1) (CT6.1) (L antigen family member 2) (LAGE-2), CTAG1A; CTAG1B, CTAG CTAG1 ESO1 LAGE2 LAGE2A; LAGE2B
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CTAG1A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.