



Klk8 Antibody

Product Code	CSB-PA720268ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q61955
Immunogen	Recombinant Mus musculus Klk8 protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-8 (mK8) (EC 3.4.21.118) (Neuropsin) (NP) (Serine protease 19), Klk8, Nrpn Prss19
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Klk8
Product Modify	Unmodified