



luxS Antibody

Product Code	CSB-PA518511ZA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O34667
Immunogen	Recombinant Bacillus subtilis luxS protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-ribosylhomocysteine lyase (EC 4.4.1.21) (AI-2 synthesis protein) (Autoinducer-2 production protein LuxS), luxS, ytiB
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	LUXS
Product Modify	Unmodified