




ASGR2 Antibody

Product Code	CSB-PA002208ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24721
Immunogen	Recombinant Mouse Asialoglycoprotein receptor 2 (Asgr2), Partial (80-301AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Asialoglycoprotein receptor 2 (ASGP-R 2) (ASGPR 2) (Hepatic lectin 2) (HL-2) (mHL-2), Asgr2, Asgr-2
Immunogen Species	Mus musculus (Mouse)
Target Names	ASGR2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ASGR2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45.9 kDa

Observed band size: 52 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.