



Three-finger toxin 3FTx-Oxy6 Antibody

Product Code	CSB-PA415963ZA01OGE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7X4T2
Immunogen	Recombinant Oxyuranus microlepidotus Three-finger toxin 3FTx-Oxy6 protein
Raised In	Rabbit
Species Reactivity	Oxyuranus microlepidotus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin 3FTx-Oxy6
Immunogen Species	Oxyuranus microlepidotus (Inland taipan)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.