



SFTPb Antibody

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA021173ZA01BO1 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P15781 |
| Immunogen | Recombinant Bos taurus SFTPb protein |
| Raised In | Rabbit |
| Species Reactivity | Bos taurus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pulmonary surfactant-associated protein B (SP-B) (6 kDa protein) (Pulmonary surfactant-associated proteolipid SPL(Phe)), SFTPb, SFTP3 |
| Immunogen Species | Bos taurus (Bovine) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | SFTPb |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |