



plcH Antibody

Product Code	CSB-PA361746ZA01EZX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06200
Immunogen	Recombinant Pseudomonas aeruginosa plcH protein
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemolytic phospholipase C (EC 3.1.4.3) (Heat-labile hemolysin) (PLC-H) (Phosphatidylcholine cholinephosphohydrolase), plcH
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG)
Lead Time	Made-to-order (12-14 weeks)
Target Names	plcH
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.