



nadB Antibody

Product Code	CSB-PA527569ZA01FHX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O57765
Immunogen	Recombinant Pyrococcus horikoshii nadB protein
Raised In	Rabbit
Species Reactivity	Pyrococcus horikoshii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-aspartate oxidase (LASPO) (EC 1.4.3.16) (Quinolinate synthase B), nadB
Immunogen Species	Pyrococcus horikoshii (strain ATCC 700860 / DSM 12428 / JCM 9974 / NBRC 100139?)
Lead Time	Made-to-order (12-14 weeks)
Target Names	nadB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.