



cut7 Antibody

Product Code	CSB-PA351343ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WP43
Immunogen	Recombinant Mycobacterium tuberculosis cut7 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable cutinase Rv1984c (EC 3.1.1.74)
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Lead Time	Made-to-order (12-14 weeks)
Target Names	cut7
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.