



Alpha-amylase inhibitor HOE-467A Antibody

Product Code	CSB-PA365519ZA01SQB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01092
Immunogen	Recombinant Streptomyces tendae Alpha-amylase inhibitor HOE-467A protein
Raised In	Rabbit
Species Reactivity	Streptomyces tendae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-amylase inhibitor HOE-467A (Tendamistat)
Immunogen Species	Streptomyces tendae
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified