



SLC4A1 Antibody

Product Code	CSB-PA021663ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02730
Immunogen	Recombinant Human Band 3 anion transport protein (SLC4A1), partial (1-403AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Band 3 anion transport protein (Anion exchange protein 1) (AE 1) (Anion exchanger 1) (Solute carrier family 4 member 1) (CD antigen CD233), SLC4A1, AE1 DI EPB3
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	SLC4A1
Product Modify	Unmodified