



Major capsid protein VP1 Antibody

Product Code	CSB-PA360953ZA01BGY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03088
Immunogen	Recombinant BK polyomavirus Major capsid protein VP1 protein
Raised In	Rabbit
Species Reactivity	BK polyomavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein VP1 (Major structural protein VP1)
Immunogen Species	BK polyomavirus (BKPyV) (Human polyomavirus 1)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.