



Mcpt2 Antibody

Product Code	CSB-PA360441ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00770
Immunogen	Recombinant Rattus norvegicus Mcpt2 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mast cell protease 2 (rMCP-2) (EC 3.4.21.-) (Group-specific protease) (Mast cell protease II) (rMCP-II), Mcpt2
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Mcpt2
Product Modify	Unmodified