



VP35 Antibody

Product Code	CSB-PA764950ZA01ZAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6V1Q9
Immunogen	Recombinant Zaire ebolavirus VP35 protein
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymerase cofactor VP35, VP35
Immunogen Species	Zaire ebolavirus (strain Kikwit-95) (ZEBOV) (Zaire Ebola virus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	VP35
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.