



fbpB Antibody

Product Code	CSB-PA314366ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WQP0
Immunogen	Recombinant Mycobacterium tuberculosis fbpB protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol acyltransferase/mycolyltransferase Ag85B (DGAT) (EC 2.3.1.122) (EC 2.3.1.20) (30 kDa extracellular protein) (Acyl-CoA:diacylglycerol acyltransferase) (Antigen 85 complex B) (85B) (Ag85B) (Extracellular alpha-antigen) (Fibronectin-binding protein B) (Fbps B), fbpB
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fbpB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.