



clpP2 Antibody

Product Code	CSB-PA528713ZA01BUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O51698
Immunogen	Recombinant Borrelia burgdorferi clpP2 protein
Raised In	Rabbit
Species Reactivity	Borrelia burgdorferi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent Clp protease proteolytic subunit 2 (EC 3.4.21.92) (Endopeptidase Clp 2), clpP2, clpP-2
Immunogen Species	Borrelia burgdorferi (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)
Lead Time	Made-to-order (12-14 weeks)
Target Names	clpP2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.