



Tick anticoagulant peptide Antibody

Product Code	CSB-PA325525ZA01OCK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17726
Immunogen	Recombinant Ornithodoros moubata Tick anticoagulant peptide protein
Raised In	Rabbit
Species Reactivity	Ornithodoros moubata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tick anticoagulant peptide (TAP)
Immunogen Species	Ornithodoros moubata (Soft tick) (Argasid tick)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.