



6GAL Antibody

Product Code	CSB-PA762323ZA01THZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76FP5
Immunogen	Recombinant Hypocrea rufa 6GAL protein
Raised In	Rabbit
Species Reactivity	Hypocrea rufa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endo-beta-1,6-galactanase (EC 3.2.1.164), 6GAL
Immunogen Species	Hypocrea rufa (Trichoderma viride)
Lead Time	Made-to-order (12-14 weeks)
Target Names	6GAL
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.