



Clotting factor C Antibody

Product Code	CSB-PA335355ZA01TAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28175
Immunogen	Recombinant Tachypleus tridentatus Clotting factor C protein
Raised In	Rabbit
Species Reactivity	Tachypleus tridentatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Limulus clotting factor C (FC) (EC 3.4.21.84) [Cleaved into: Limulus clotting factor C heavy chain; Limulus clotting factor C light chain; Limulus clotting factor C chain A; Limulus clotting factor C chain B]
Immunogen Species	Tachypleus tridentatus (Japanese horseshoe crab)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.