



glcB Antibody

Product Code	CSB-PA403987ZA01MON
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5U3K4
Immunogen	Recombinant Mycobacterium tuberculosis glcB protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Malate synthase G (EC 2.3.3.9), glcB
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25177 / H37Ra)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	glcB
Product Modify	Unmodified