



ligA Antibody

Product Code	CSB-PA669648ZA01CBAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q47YI0
Immunogen	Recombinant Colwellia psychrerythraea ligA protein
Raised In	Rabbit
Species Reactivity	Colwellia psychrerythraea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA ligase (EC 6.5.1.2) (Polydeoxyribonucleotide synthase [NAD(+)]), ligA
Immunogen Species	Colwellia psychrerythraea (strain 34H / ATCC BAA-681) (Vibrio psychroerythus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ligA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.