



# S100A9 Antibody

<b>Product Code</b>	CSB-PA020642ZA01BO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P28783
<b>Immunogen</b>	Recombinant Bos taurus S100A9 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bos taurus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein S100-A9 (BEE22) (Calgranulin-B) (Neutrophil cytosolic 23 kDa protein) (p23) (S100 calcium-binding protein A9), S100A9
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	S100A9
<b>Product Modify</b>	Unmodified