



eltA Antibody

Product Code	CSB-PA356922ZA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06717
Immunogen	Recombinant Escherichia coli Heat-labile enterotoxin A chain (eltA) (19-258AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-labile enterotoxin A chain (LT-A, porcine) (LTP-A), eltA, ltpA
Immunogen Species	Escherichia coli
Lead Time	Made-to-order
Gene Names	eltA
Product Modify	Unmodified