



Mafk Antibody

Product Code	CSB-PA713990ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q61827
Immunogen	Recombinant Mouse Transcription factor MafK (Mafk) (1-156AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor MafK (Erythroid transcription factor NF-E2 p18 subunit), Mafk, Nfe2u
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Gene Names	Mafk
Product Modify	Unmodified