



ABCC1 Antibody

Product Code	CSB-PA001056ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33527
Immunogen	Recombinant Human Multidrug resistance-associated protein 1 (ABCC1), partial (1248-1531AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Multidrug resistance-associated protein 1 (ATP-binding cassette sub-family C member 1) (Leukotriene C(4) transporter) (LTC4 transporter), ABCC1, MRP MRP1
Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	ABCC1
Product Modify	Unmodified