



TK Antibody

Product Code	CSB-PA023577ZA01HWY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03176
Immunogen	Recombinant Human herpesvirus 1 TK protein
Raised In	Rabbit
Species Reactivity	Human herpesvirus 1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymidine kinase (EC 2.7.1.21), TK
Immunogen Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	TK
Product Modify	Unmodified