



E7 Antibody

Product Code	CSB-PA339740ZA01HNY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36831
Immunogen	Recombinant Human papillomavirus type 52 E7 protein
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 52
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E7, E7
Immunogen Species	Human papillomavirus type 52
Lead Time	Made-to-order (12-14 weeks)
Target Names	E7
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.