



E7 Antibody

Product Code	CSB-PA365855ZA01HML
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03129
Immunogen	Recombinant Human papillomavirus type 16 Protein E7 (E7) (1-98AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 16
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E7, E7
Immunogen Species	Human papillomavirus type 16
Lead Time	Made-to-order
Gene Names	E7
Product Modify	Unmodified