



F Antibody

Product Code	CSB-PA356041ZA01HPO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03420
Immunogen	Recombinant Human respiratory syncytial virus A F protein
Raised In	Rabbit
Species Reactivity	Human respiratory syncytial virus A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fusion glycoprotein F0 (Protein F) [Cleaved into: Fusion glycoprotein F2'; Interchain peptide; Fusion glycoprotein F2; Fusion glycoprotein F1], F
Immunogen Species	Human respiratory syncytial virus A (strain A2)
Lead Time	Made-to-order (12-14 weeks)
Target Names	F
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.