



BH0637 Antibody

Product Code	CSB-PA864561ZA01BQS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9KF49
Immunogen	Recombinant Bacillus halodurans BH0637 protein
Raised In	Rabbit
Species Reactivity	Bacillus halodurans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative adenine deaminase BH0637 (Adenase) (Adenine aminase) (EC 3.5.4.2)
Immunogen Species	Bacillus halodurans (strain ATCC BAA-125 / DSM 18197 / FERM 7344 / JCM 9153 / C-125)
Lead Time	Made-to-order (12-14 weeks)
Target Names	BH0637
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.