



lytA Antibody

Product Code	CSB-PA326479ZA01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24556
Immunogen	Recombinant Staphylococcus aureus lytA protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Autolysin (EC 3.5.1.28) (N-acetylmuramoyl-L-alanine amidase), lytA
Immunogen Species	Staphylococcus aureus
Lead Time	Made-to-order (12-14 weeks)
Target Names	lytA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.