



SAA1 Antibody

Product Code	CSB-PA020656ZA01CA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19707
Immunogen	Recombinant Felis catus SAA1 protein
Raised In	Rabbit
Species Reactivity	Felis catus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum amyloid A protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA)], SAA1
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	SAA1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.