



hph Antibody

Product Code	CSB-PA365448ZA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00557
Immunogen	Recombinant Escherichia coli Hygromycin-B 4-O-kinase (hph) (1-341AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hygromycin-B 4-O-kinase (EC 2.7.1.163) (APH(4)) (Hygromycin B phosphotransferase) (Hygromycin-B kinase), hph
Immunogen Species	Escherichia coli
Lead Time	Made-to-order
Gene Names	hph
Product Modify	Unmodified