




# F17a-A Antibody

<b>Product Code</b>	CSB-PA318126ZA01ENL
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P11312
<b>Immunogen</b>	Recombinant Escherichia coli F17 fimbrial protein (F17a-A) (22-180AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Tested Applications</b>	ELISA, WB
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein A/G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	F17 fimbrial protein (F17 fimbrial adhesive antigen) (F17 pilin), F17a-A
<b>Immunogen Species</b>	Escherichia coli
<b>Target Names</b>	F17a-A

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: F17a-A antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 32.4 kDa

Observed band size: 38 kDa

<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.