




Trypsin/alpha-amylase inhibitor CMX1/CMX3 Antibody

Product Code	CSB-PA675990ZA01TQN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q43723
Immunogen	Recombinant Triticum aestivum Trypsin/alpha-amylase inhibitor CMX1/CMX3 protein (25-121aa)
Raised In	Rabbit
Species Reactivity	Triticum aestivum
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Triticum aestivum
Target Names	Trypsin/alpha-amylase inhibitor CMX1/CMX3

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: Trypsin/alpha-amylase inhibitor

CMX1/CMX3 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 27.5 kDa

Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.