



Vwf Antibody

Product Code	CSB-PA025960ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8CIZ8
Immunogen	Recombinant Mus musculus Vwf protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	von Willebrand factor (vWF) [Cleaved into: von Willebrand antigen 2 (von Willebrand antigen II)], Vwf
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Vwf
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.