



gB Antibody

Product Code	CSB-PA309330ZA01RKC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P89053
Immunogen	Recombinant Rhesus cytomegalovirus gB protein
Raised In	Rabbit
Species Reactivity	Rhesus cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein B (gB), gB
Immunogen Species	Rhesus cytomegalovirus (strain 68-1) (RhCMV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	gB
Product Modify	Unmodified