



folP Antibody

Product Code	CSB-PA682423ZA01SLF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59919
Immunogen	Recombinant Staphylococcus haemolyticus folP protein
Raised In	Rabbit
Species Reactivity	Staphylococcus haemolyticus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropteroate synthase (DHPS) (EC 2.5.1.15) (Dihydropteroate pyrophosphorylase), folP
Immunogen Species	Staphylococcus haemolyticus
Lead Time	Made-to-order (12-14 weeks)
Target Names	folP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.