



ACT1 Antibody

Product Code	CSB-PA320814ZA01AAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10982
Immunogen	Recombinant Absidia glauca ACT1 protein
Raised In	Rabbit
Species Reactivity	Absidia glauca
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-1 (Fragment), ACT1
Immunogen Species	Absidia glauca (Pin mould)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ACT1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.