



ospA Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA363424ZA01BIO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0A3N6 |
| Immunogen | Recombinant Borreliella burgdorferi ospA protein |
| Raised In | Rabbit |
| Species Reactivity | Borreliella burgdorferi |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Outer surface protein A, ospA |
| Immunogen Species | Borreliella burgdorferi (Lyme disease spirochete) (Borrelia burgdorferi) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | ospA |
| Product Modify | Unmodified |