



P39 Antibody

Product Code	CSB-PA327448ZA01LRP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35840
Immunogen	Recombinant Lymantria dispar multicapsid nuclear polyhedrosis virus P39 protein
Raised In	Rabbit
Species Reactivity	Lymantria dispar multicapsid nuclear polyhedrosis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major capsid protein, P39
Immunogen Species	Lymantria dispar multicapsid nuclear polyhedrosis virus (LdMNPV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	P39
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.