



p37 Antibody

Product Code	CSB-PA675437ZA01MLN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q49410
Immunogen	Recombinant Mycoplasma genitalium p37 protein
Raised In	Rabbit
Species Reactivity	Mycoplasma genitalium
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity transport system protein p37, p37
Immunogen Species	Mycoplasma genitalium (strain ATCC 33530 / G-37 / NCTC 10195)
Lead Time	Made-to-order (12-14 weeks)
Target Names	p37
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.