




Ang4 Antibody

Product Code	CSB-PA661010ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3TMQ6
Immunogen	Recombinant Mus musculus (Mouse) Angiogenin-4 (25-144aa)
Species Reactivity	Mus musculus (Mouse)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	Ang4

Image



Western Blot
Positive WB detected in Recombinant protein
All lanes: Ang4 antibody at 1:2000
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 31.4 kDa
Observed band size: 33 kDa