



GHR Antibody

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA009411ZA01MOW |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P79194 |
| Immunogen | Recombinant Macaca mulatta GHR protein |
| Raised In | Rabbit |
| Species Reactivity | Macaca mulatta |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Growth hormone receptor (GH receptor) (Somatotropin receptor) [Cleaved into: Growth hormone-binding protein (GH-binding protein) (GHBP) (Serum-binding protein)], GHR |
| Immunogen Species | Macaca mulatta (Rhesus macaque) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | GHR |
| Product Modify | Unmodified |